# ANTI-AGE

MAGAZINE

## Successful Aging

25 experts vous conseillent  
Advice from 25 experts

**Votre nutrition  
anti-âge**  
Your anti-aging  
nutrition

**Anti Taches  
et Photo-  
vieillesse**  
Fighting dark  
marks and  
photoaging

**Découvrir  
la cryolipolyse**  
Everything  
you need to know  
about cryolipolysis

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## NEAUVIA STIMULATE: A NEW POLYMERIC MATRIX ENRICHED WITH CAHA

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
Neauvia Stimulate is a PEG-crosslinked hydrogel with 1% calcium hydroxyapatite microspheres trapped in the mesh of the polymeric matrix (Fig. 1), which prevents them from migrating or conglomerating after implantation.

Neauvia Hydrogel Matrix (NHM) is highly biocompatible, so there is only a mild inflammatory response and slow enzymatic degradation.

The reduced degradation via enzymatic mechanisms and mechanical stress promotes an immediate and longer lasting corrective effect. In addition, the modulated interaction of the lysosomes in the cell's phagocytic system with the small hydroxyapatite calcium microspheres (10-12 micron) helps with stimulating fibroblasts to produce new extracellular matrix and collagen in particular. The newly formed collagen, made up of fine, type III collagen fibrils, integrates harmoniously with the molecular elements of the extracellular matrix ground substance, with an effective additional supporting, filling and revitalizing action.

In vitro studies have documented the capacity of the new PEG-crosslinked filler to induce metabolic activation in fibroblasts, even at low concentrations of CaHA.

Fig. 2

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*Histological sections stained by Picrosirius Red and observed by circularly polarized light microscopy. In both figures, red/orange areas represent old collagen fibres and green/yellow areas represent newly formed fibres. In the skin sample following Stimulate treatment, it is evident the high presence of green/yellow collagen fibres, representing newly formed fibres (collagen III).*

Avant / Before

Après / After

The reason behind this action probably lies in NHM's high level of biointegration, which allows to create three-dimensional interstitial spaces that are fully integrated with the cellular and fibrillar structures in connective tissues.

This complete biointegration, with reduced inflammatory activity, is confirmed in vitro and in vivo. It allows Neauvia Stimulate to combine its high filling capacity (thanks to the hydrogel's excellent rheological properties) with a documented stimulating action on the neocollagenesis process (Fig. 2).

The clinical indications of Neauvia Stimulate include the correction of deep flaws and, additionally, a further interaction through a biostimulating (non mechanical) pathway that involves the more superficial elements, regenerating also the fibrillar components of the dermis.

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Degree in General Medicine and Surgery from the University of Pavia. Postgraduate Degree in Dermatology and Venereology. Teacher of Aesthetic Medicine for the Postgraduate Course in Aesthetic medicine at the University of Pavia, directed by Prof. G. Robutti. Dr. Zerbinati works as dermatologist in his private practices - CMP and since November 2000 teaches as full Professor of Dermatology and Venereology at the University of Insubria, Varese (Italy). Scientific Coordinator of the International Academy of Cosmetic Gynecology IACG.

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(toxine botulique, fillers, laser, lumière pulsée, etc.), dont elle optimise les résultats. mesohyal® is the widest range of class III medical devices with CE mark that improves the appearance of the skin immediately, naturally and providing long lasting effect. The range is ideal for combining with and optimising the results of other medical procedures (botulinum toxin, fillers, laser, pulsed light, etc.).

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pour le rajeunissement des mains. Neuvia offers a biorevitalizer (Hyaluronic acid without cross-linking PEG) HYDRO DELUXE: a soft treatment for face and body with hyaluronic acid, which in 4/5 sessions, hydrates and re-densifies tissues by bio-stimulating collagen (contains 0,1% of calcium hydroxyapatite). Ideal for the biorevitalization of delicate areas around the eyes and neck and for hand rejuvenation.

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